

## Product datasheet for **GA109893**

### Human Lipolysis Stimulated Lipoprotein Receptor (LSR) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	LSR CRISPRa kit - CRISPR gene activation of human lipolysis stimulated lipoprotein receptor
Symbol:	LSR
UniProt ID:	<a href="#">Q86X29</a>
Synonyms:	ILDR3; LISCH7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)