

Product datasheet for **GA109908**

Human DC2L1 (DYNC2LI1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DYNC2LI1 CRISPRa kit - CRISPR gene activation of human dynein cytoplasmic 2 light intermediate chain 1
Symbol:	DYNC2LI1
UniProt ID:	Q8TCX1
Synonyms:	CGI-60; D2LIC; DKFZp564A033; dynein, cytoplasmic 2, light intermediate chain 1; dynein 2 light intermediate chain; LIC3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)