

## Product datasheet for **GA109915**

### Human SF3B14 (SF3B6) activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | SF3B6 CRISPRa kit - CRISPR gene activation of human splicing factor 3b subunit 6 |
| Symbol:       | SF3B6  |
| UniProt ID:   | <a href="#">Q9Y3B4</a>   |
| Synonyms:     | CGI-110; HSPC175; Ht006; P14; SAP14; SAP14a; SF3B14; SF3B14a                     |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)          |

#### Product images:



[View online >](#)