

Product datasheet for **GA109936**

Human TMEM66 (SARAF) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SARAF CRISPRa kit - CRISPR gene activation of human store-operated calcium entry associated regulatory factor
Symbol:	SARAF
Components:	GA109936G1, TMEM66 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109936G2, TMEM66 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109936G3, TMEM66 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96BY9</u>
Synonyms:	FOAP-7; HSPC035; TMEM66; XTP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

