

Product datasheet for **GA109937**

Human TPPP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TPPP3 CRISPRa kit - CRISPR gene activation of human tubulin polymerization promoting protein family member 3
Symbol:	TPPP3
Components:	GA109937G1, TPPP3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109937G2, TPPP3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109937G3, TPPP3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9BW30
Synonyms:	CGI-38; p20; p25gamma; TPPP/p20
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)