

## Product datasheet for **GA109944**

### Human CPSF73 (CPSF3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	CPSF3 CRISPRa kit - CRISPR gene activation of human cleavage and polyadenylation specific factor 3
Symbol:	CPSF3
UniProt ID:	<a href="#">Q9UKF6</a>
Synonyms:	CPSF-73; CPSF73
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)