

Product datasheet for **GA109946**

Human HECA activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HECA CRISPRa kit – CRISPR gene activation of human hdc homolog, cell cycle regulator
Symbol:	HECA
Components:	GA109946G1, HECA gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109946G2, HECA gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109946G3, HECA gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UBI9</u>
Synonyms:	dJ225E12.1; HDC; HDCL; HHDC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

