

Product datasheet for **GA109950**

Human PADI3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PADI3 CRISPRa kit - CRISPR gene activation of human peptidyl arginine deiminase 3
Symbol:	PADI3
UniProt ID:	Q9ULW8
Synonyms:	PAD3; PDI3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)