


Product datasheet for **GA109955**

Human ZNF44 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF44 CRISPRa kit - CRISPR gene activation of human zinc finger protein 44
Symbol:	ZNF44
Components:	GA109955G1, ZNF44 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109955G2, ZNF44 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109955G3, ZNF44 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P15621</u>
Synonyms:	GIOT-2; KOX7; ZNF; ZNF55; ZNF58; ZNF504
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)