

Product datasheet for **GA109959**

Human RAP80 (UIMC1) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | UIMC1 CRISPRa kit - CRISPR gene activation of human ubiquitin interaction motif containing 1 |
| Symbol: | UIMC1 |
| UniProt ID: | <u>Q96RL1</u> |
| Synonyms: | RAP80; X2HRIP110 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)