

Product datasheet for **GA109979**

Human SYT17 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SYT17 CRISPRa kit - CRISPR gene activation of human synaptotagmin 17
Symbol:	SYT17
Components:	GA109979G1, SYT17 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109979G2, SYT17 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109979G3, SYT17 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BSW7</u>
Synonyms:	sytXVII
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)