

Product datasheet for **GA109985**

Human TM7SF3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TM7SF3 CRISPRa kit - CRISPR gene activation of human transmembrane 7 superfamily member 3
Symbol:	TM7SF3
Components:	GA109985G1, TM7SF3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109985G2, TM7SF3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109985G3, TM7SF3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NS93</u>
Synonyms:	2010003B14Rik; NARP1; OTTMUSP00000031284; transmembrane 7 superfamily member 3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)