

Product datasheet for **GA109995**

Human PHAX activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PHAX CRISPRa kit - CRISPR gene activation of human phosphorylated adaptor for RNA export
Symbol:	PHAX
Components:	GA109995G1, PHAX gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109995G2, PHAX gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109995G3, PHAX gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9H814
Synonyms:	RNUXA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

