

Product datasheet for **GA110006**

Human **TM6SF1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TM6SF1 CRISPRa kit - CRISPR gene activation of human transmembrane 6 superfamily member 1
Symbol:	TM6SF1
Components:	GA110006G1, TM6SF1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110006G2, TM6SF1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110006G3, TM6SF1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9BZW5
Synonyms:	transmembrane 6 superfamily member 1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)