

Product datasheet for **GA110011**

Human SHC3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SHC3 CRISPRa kit - CRISPR gene activation of human SHC adaptor protein 3
Symbol:	SHC3
UniProt ID:	Q92529
Synonyms:	N-Shc; NSHC; RAI; SHCC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)