

## Product datasheet for **GA110013**

### Human TPCN1 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TPCN1 CRISPRa kit - CRISPR gene activation of human two pore segment channel 1
Symbol:	TPCN1
UniProt ID:	<u><a href="#">Q9ULQ1</a></u>
Synonyms:	TPC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)