

Product datasheet for **GA110019**

Human AMPK gamma 3 (PRKAG3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRKAG3 CRISPRa kit - CRISPR gene activation of human protein kinase AMP-activated non-catalytic subunit gamma 3
Symbol:	PRKAG3
UniProt ID:	Q9UGI9
Synonyms:	AMPKG3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)