

Product datasheet for **GA110020**

Human PTOV1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PTOV1 CRISPRa kit - CRISPR gene activation of human PTOV1 extended AT-hook containing adaptor protein
Symbol:	PTOV1
Components:	GA110020G1, PTOV1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110020G2, PTOV1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110020G3, PTOV1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86YD1</u>
Synonyms:	ACID2; PTOV-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)