

Product datasheet for **GA110027**

Human GPR86 (P2RY13) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2RY13 CRISPRa kit - CRISPR gene activation of human purinergic receptor P2Y13
Symbol:	P2RY13
UniProt ID:	Q9BPV8
Synonyms:	FKSG77; G-protein coupled receptor 94; GPCR1; GPCR1, GPR86, GPR94, P2Y13, SP174, FKSG77; GPR86; GPR94; G protein-coupled receptor 86; P2Y13; P2Y purinoceptor 13; purinergic receptor P2Y, G-protein coupled, 13; SP174
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)