


Product datasheet for **GA110043**

Human PPIL3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPIL3 CRISPRa kit - CRISPR gene activation of human peptidylprolyl isomerase like 3
Symbol:	PPIL3
Components:	GA110043G1, PPIL3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110043G2, PPIL3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110043G3, PPIL3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H2H8</u>
Synonyms:	CYPJ
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)