

Product datasheet for **GA110050**

Human glycerol 3 phosphate permease (SLC37A1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC37A1 CRISPRa kit - CRISPR gene activation of human solute carrier family 37 member 1
Symbol:	SLC37A1
UniProt ID:	P57057
Synonyms:	G3PP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)