

Product datasheet for **GA110052**

Human PCBP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCBP3 CRISPRa kit - CRISPR gene activation of human poly(rC) binding protein 3
Symbol:	PCBP3
Components:	GA110052G1, PCBP3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110052G2, PCBP3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110052G3, PCBP3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	P57721
Synonyms:	ALPHA-CP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

