


Product datasheet for **GA110057**

Human TSPEAR activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSPEAR CRISPRa kit - CRISPR gene activation of human thrombospondin type laminin G domain and EAR repeats
Symbol:	TSPEAR
Components:	GA110057G1, TSPEAR gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110057G2, TSPEAR gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110057G3, TSPEAR gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8WU66</u>
Synonyms:	C21orf29; DFNB98; TSP-EAR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)