

Product datasheet for **GA110074**

Human KCNK10 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNK10 CRISPRa kit - CRISPR gene activation of human potassium two pore domain channel subfamily K member 10
Symbol:	KCNK10
UniProt ID:	P57789
Synonyms:	K2p10.1; PPP1R97; TREK-2; TREK2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)