

Product datasheet for **GA110079**

Human SREB3 (GPR173) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GPR173 CRISPRa kit - CRISPR gene activation of human G protein-coupled receptor 173
Symbol:	GPR173
UniProt ID:	Q9NS66
Synonyms:	SREB3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)