

Product datasheet for **GA110085**

Human UNC93A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	UNC93A CRISPRa kit - CRISPR gene activation of human unc-93 homolog A
Symbol:	UNC93A
Components:	GA110085G1, UNC93A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110085G2, UNC93A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110085G3, UNC93A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q86WB7
Synonyms:	dj366N23.1; dj366N23.2; Unc-93A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)