

Product datasheet for **GA110089**

Human RAP1 (TERF2IP) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TERF2IP CRISPRa kit - CRISPR gene activation of human TERF2 interacting protein
Symbol:	TERF2IP
UniProt ID:	<u>Q9NYB0</u>
Synonyms:	DRIP5; RAP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)