


Product datasheet for **GA110090**

Human SLC38A2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC38A2 CRISPRa kit - CRISPR gene activation of human solute carrier family 38 member 2
Symbol:	SLC38A2
Components:	GA110090G1, SLC38A2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110090G2, SLC38A2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110090G3, SLC38A2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96QD8</u>
Synonyms:	ATA2; PRO1068; SAT2; SNAT2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)