


Product datasheet for **GA110102**

Human SASH3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SASH3 CRISPRa kit - CRISPR gene activation of human SAM and SH3 domain containing 3
Symbol:	SASH3
Components:	GA110102G1, SASH3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110102G2, SASH3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110102G3, SASH3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O75995</u>
Synonyms:	753P9; CXorf9; HAC2; SH3D6C; SLY
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)