

Product datasheet for **GA110109**

Human TAF7L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAF7L CRISPRa kit - CRISPR gene activation of human TATA-box binding protein associated factor 7 like
Symbol:	TAF7L
UniProt ID:	Q5H9L4
Synonyms:	CT40; TAF2Q
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)