

Product datasheet for **GA110113**

Human **FAM190B (CCSER2)** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	CCSER2 CRISPRa kit - CRISPR gene activation of human coiled-coil serine rich protein 2
Symbol:	CCSER2
Components:	GA110113G1, FAM190B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110113G2, FAM190B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110113G3, FAM190B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H7UI</u>
Synonyms:	bA486O22.1; FAM190B; Gcap14; KIAA1128
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

