


## Product datasheet for **GA110117**

### Human SPIN2A activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | SPIN2A CRISPRa kit - CRISPR gene activation of human spindlin family member 2A                                                                                                                                                                                                                                |
| Symbol:       | SPIN2A                                                                                                                                                                                                                                                                                                        |
| Components:   | GA110117G1, SPIN2A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA110117G2, SPIN2A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA110117G3, SPIN2A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q99865</a>                                                                                                                                                                                                                                                                                        |
| Synonyms:     | dj323P24.1; DXF34; SPIN2; TDRD25                                                                                                                                                                                                                                                                              |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)