

Product datasheet for **GA110127**

Human PEPP2 (PLEKHA5) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLEKHA5 CRISPRa kit - CRISPR gene activation of human pleckstrin homology domain containing A5
Symbol:	PLEKHA5
Components:	GA110127G1, PEPP2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110127G2, PEPP2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110127G3, PEPP2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9HAU0</u>
Synonyms:	PEPP-2; PEPP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)