

Product datasheet for **GA110136**

Human TMX3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMX3 CRISPRa kit - CRISPR gene activation of human thioredoxin related transmembrane protein 3
Symbol:	TMX3
Components:	GA110136G1, TMX3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110136G2, TMX3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110136G3, TMX3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96JJ7
Synonyms:	PDIA13; TXNDC10
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)