

Product datasheet for **GA110150**

Human DDX4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DDX4 CRISPRa kit - CRISPR gene activation of human DEAD-box helicase 4
Symbol:	DDX4
UniProt ID:	<u>Q9NQI0</u>
Synonyms:	DEAD (Asp-Glu-Ala-Asp) box polypeptide 4; DEAD/H (Asp-Glu-Ala-Asp/His) box polypeptide 4; MGC111074; OTTHUMP00000122546; VASA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)