

Product datasheet for **GA110155**

Human WDR44 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR44 CRISPRa kit - CRISPR gene activation of human WD repeat domain 44
Symbol:	WDR44
Components:	GA110155G1, WDR44 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110155G2, WDR44 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110155G3, WDR44 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5JSH3
Synonyms:	RAB11BP; RPH11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)