

Product datasheet for **GA110159**

Human USP53 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	USP53 CRISPRa kit - CRISPR gene activation of human ubiquitin specific peptidase 53
Symbol:	USP53
Components:	GA110159G1, USP53 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110159G2, USP53 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110159G3, USP53 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q70EK8</u>
Synonyms:	DKFZp781E1417
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)