

Product datasheet for **GA110160**

Human MRPL50 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MRPL50 CRISPRa kit - CRISPR gene activation of human mitochondrial ribosomal protein L50
Symbol:	MRPL50
Components:	GA110160G1, MRPL50 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110160G2, MRPL50 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110160G3, MRPL50 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N5N7</u>
Synonyms:	MRP-L50
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)