

Product datasheet for **GA110163**

Human **FAM35A** activation kit by **CRISPRa**

Product data:

| | |
|---------------|--------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits |
| Description: | SHLD2 CRISPRa kit - CRISPR gene activation of human shieldin complex subunit 2 |
| Symbol: | SHLD2 |
| UniProt ID: | <u>Q86V20</u> |
| Synonyms: | bA163M19.1; FAM35A1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)