

Product datasheet for **GA110194**

Human L1TD1 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | L1TD1 CRISPRa kit - CRISPR gene activation of human LINE1 type transposase domain containing 1 |
| Symbol: | L1TD1 |
| UniProt ID: | Q5T7N2 |
| Synonyms: | ECAT11 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)