

Product datasheet for **GA110203**

Human PAF1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PAF1 CRISPRa kit - CRISPR gene activation of human PAF1 homolog, Paf1/RNA polymerase II complex component
Symbol:	PAF1
UniProt ID:	<u>Q8N7H5</u>
Synonyms:	F23149_1; PD2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)