

## Product datasheet for **GA110215**

### Human LRRN3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	LRRN3 CRISPRa kit - CRISPR gene activation of human leucine rich repeat neuronal 3
Symbol:	LRRN3
UniProt ID:	<a href="#">Q9H3W5</a>
Synonyms:	FIGLER5; NLRR-3; NLRR3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)