


Product datasheet for **GA110222**

Human **RRN3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	RRN3 CRISPRa kit - CRISPR gene activation of human RRN3 homolog, RNA polymerase I transcription factor
Symbol:	RRN3
Components:	GA110222G1, RRN3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110222G2, RRN3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110222G3, RRN3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NYV6</u>
Synonyms:	A-270G1.2; TIFIA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)