


Product datasheet for **GA110233**

Human RAB39 (RAB39A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RAB39A CRISPRa kit - CRISPR gene activation of human RAB39A, member RAS oncogene family
Symbol:	RAB39A
Components:	GA110233G1, RAB39 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110233G2, RAB39 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110233G3, RAB39 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q14964
Synonyms:	RAB39
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)