

Product datasheet for **GA110237**

Human LEPROT activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LEPROT CRISPRa kit - CRISPR gene activation of human leptin receptor overlapping transcript
Symbol:	LEPROT
UniProt ID:	O15243
Synonyms:	LEPR; OB-RGRP; OBRGRP; VPS55
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)