

Product datasheet for **GA110244**

Human IL17RD activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IL17RD CRISPRa kit - CRISPR gene activation of human interleukin 17 receptor D
Symbol:	IL17RD
UniProt ID:	<u>Q8NFM7</u>
Synonyms:	HH18; IL-17RD; IL17RLM; SEF
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)