

Product datasheet for **GA110245**

Human FA20A (FAM20A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FAM20A CRISPRa kit - CRISPR gene activation of human FAM20A golgi associated secretory pathway pseudokinase
Symbol:	FAM20A
Components:	GA110245G1, FA20A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110245G2, FA20A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110245G3, FA20A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96MK3
Synonyms:	AI1G; AIGFS; FP2747
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)