

## Product datasheet for **GA110247**

### Human Proprotein convertase PC4 (PCSK4) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PCSK4 CRISPRa kit - CRISPR gene activation of human proprotein convertase subtilisin/kexin type 4
Symbol:	PCSK4
UniProt ID:	<a href="#">Q6UW60</a>
Synonyms:	PC4; SPC5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)