

Product datasheet for **GA110257**

Human RNF111 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RNF111 CRISPRa kit - CRISPR gene activation of human ring finger protein 111
Symbol:	RNF111
Components:	GA110257G1, RNF111 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110257G2, RNF111 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110257G3, RNF111 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZNA4</u>
Synonyms:	ARK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)