

Product datasheet for **GA110264**

Human TRPM4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TRPM4 CRISPRa kit - CRISPR gene activation of human transient receptor potential cation channel subfamily M member 4
Symbol:	TRPM4
UniProt ID:	Q8TD43
Synonyms:	hTRPM4; LTrpC4; PFHB1B; TRPM4B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)